

Prakriti Track

Biomedical Waste Management System — Trial Guide

This guide provides a complete, step-by-step walkthrough to test the end-to-end functionality of the Prakriti Track platform. The system is divided into three key roles: **Admin**, **Driver**, and **Plant Manager**.

1. Initial Setup (Admin)

To be performed on a desktop or laptop.

1. Navigate to the application URL in your web browser.
 2. Log in using the **Plant Head** credentials (`admin@prakrititrack.in` / `admin123`).
 3. **Register a Healthcare Facility (HCF):**
 - Go to **Management > HCF Registry**.
 - Click + **Add HCF** and create a test hospital (ensure you assign it an HCF Code, e.g., `HCF001`).
 4. **Register a Vehicle & Driver:**
 - Go to **Management > Vehicles**.
 - Click + **Add Vehicle** (e.g., `JH-10-AB-1234`) and assign `Driver` from the dropdown.
-

2. Generate Waste Bags (Admin)

Simulating the hospital handing over bags.

1. Go to **Bag Tracker** from the sidebar.
 2. Click + **Create Bags**.
 3. Select your test hospital from the dropdown.
 4. Select a waste category (e.g., **Yellow**) and enter a quantity (e.g., `3`).
 5. Click **Generate Bags**.
 6. The system will instantly generate unique QR codes for each bag. Click **Print All Labels** to view the labels.
Note: For the trial, you can print these labels, or simply leave them open on your screen to scan with a mobile phone.
-

3. Collection Workflow (Driver)

Simulating the driver picking up the waste. Perform these steps on a mobile phone or split-screen browser.

1. Log in using the **Driver** credentials (`driver1@prakrititrack.in` / `driver123`).
 2. **Start Route:**
 - On the Route Overview page, select your vehicle and click **Start Today's Route**.
 3. **Check-In:**
 - Click **GPS Check-in** to log your arrival coordinates at the hospital.
 4. **Scan Bags:**
 - Click **Scan Bag** to open the built-in camera.
 - Point the camera at the 3 QR codes you generated in Step 2.
 - After scanning a code, click **Confirm Collection**.
 5. **Weigh Bags (Optional):**
 - Go to the **Weigh** tab and manually enter the weight for the collected bags.
-

4. Plant Gate Receipt (Plant Manager)

Simulating the truck arriving at the treatment plant.

1. Log in using the **Manager** credentials (`manager@prakrititrack.in` / `manager123`).
 2. Go to **Plant > Gate Scan**.
 3. Click **Open Camera Scanner** and scan the exact same 3 QR codes to simulate unloading the truck at the facility.
 4. Once all bags appear on the list, click **Confirm Receipt**.
-

5. Reconciliation (Plant Manager)

Ensuring zero waste divergence during transit.

1. Go to **Plant > Reconciliation**.
 2. Select the active route driven by your driver.
 3. The system will compare the bags marked "Collected" against the bags marked "Received". If successfully completed, the difference will show as **0**.
-

6. Treatment & Certification (Plant Manager)

Simulating the final disposal process.

1. **Create a Batch:**
 - Go to **Plant > Batches**.
 - Select the 3 received bags from the table and click + **Create Batch**.
 - Assign a treatment type (e.g., Autoclave) and save.
 2. **Mark as Treated:**
 - Go to **Plant > Treatment & Certs**.
 - Find your newly created batch and click **View**.
 - Click **Mark as Treated**.
 3. **Generate Certificate:**
 - Click **Print Certificate** to instantly generate the official Disposal Certificate for the hospital.
-

End of Trial Workflow.